

GI- LUTE

Glass Ionomer Luting Cement

TECHNICAL & SAFETY INSTRUCTIONS FOR USE



Product Description: Glass Ionomer Luting Cement Type I is based on radiopaque fine glass powder and organic polymer liquid used as luting material with excellent tooth adhesion properties having high fluoride release and flow ability

Application:

- 1.Cementation of metal-based inlays, inlays crown and bridges.
- 2.Cementation of high strength (Zirconium based) all ceramic crowns and bridges.
- 3.Cementation of stainless steel crowns
- 4.Attachment of metal orthodontic bands

Composition: Powder - Calcium-alumina-fluorosilicate glass and Liquid - Polyacrylic acid.

Methods of use: Shake the powder bottle well every time before use. Dispense one scoop of powder and two drops of liquid on to the mixing plate (The standard powder to liquid ratio is 1.5g:1g). Use plastic spatula to divide the powder into 2 equal parts and mix the first portion with all of the liquid for 10 seconds. Mix the remaining portion with the paste formed for 20 seconds. The total mixing time should be 30 seconds. Transfer cement to the prepared tooth surface using suitable placement instrument and finish the cementation. Allow 8-10 minutes for proper Setting.

Storage: Store at room temperature between 10- 35°C . Protect from direct sunlight. Close the container carefully after use. Keep away from moisture. Shelf life not affected after opening, if stored in a cool & dry place.

Disposable conditions: To be disposed off based on the prevailing rules of the country in which it is used.


Disclaimer: The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instruction of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.














Warning: Avoid prolonged or repeated contact with skin (allergic contact dermatitis) and eyes. Avoid prolonged inhalation. Do not take it internally. Contact with the skin may cause irritation, burns of hypersensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

Presentation : 1 x 15 gm, 1 x 10 ml, RefNo. RD11-03, 1 x Mixing Pad, 1 x Spoon.

Shelf Life: 3 years from the date of Manufacturing.

NOTE: TO BE USED BY DENTAL PROFESSIONALS ONLY.

 Manufactured By:
Golden Bird International
Plot No. A-4 Vishwakarma
Industrial Estate, Village Panchli
Khurd Meerut-250002 India, info@indiagbi.com
www.dentakris.com, GOLDEN-Meeru-UP/M/MD/005962
®-Registered Trademark Only

 Symbol for "BATCH CODE"  Symbol for "USE BY"  Symbol for "NON STERILE"  Symbol for "KEEP DRY"  Symbol for "CAUTION"
 SYMBOL FOR EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE  Symbol for "KEEP AWAY FROM SUNLIGHT"
 Symbol for "KEEP OUT OF THE REACH OF CHILDREN"  Symbol for "MANUFACTURER NAME & ADDRESS"  Symbol for "CONSULT INSTRUCTION FOR USE"  Symbol for "CATALOGUE NUMBER"  Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"
 Symbol for "TEMPERATURE LIMITATION"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.