

# Insta Relief

## Pain Relief Oil

### TECHNICAL & SAFETY INSTRUCTIONS FOR USE



**INTENDED USE:** Insta Relief has long been used as a natural remedy for toothaches, and its medicinal properties have been used to treat a variety of oral ailments. In recent years, it has gained increased attention as a potential solution for reducing dental pain. explore how it works to reduce dental pain and provide tips for using this natural remedy for toothaches and other oral pain. By understanding how it can help reduce dental pain, you can better decide if it is the right choice for you.

**Composition:** It is an aromatic essential oil obtained from Insta relief have a long history of use in traditional medicine, both in extract and essential oil form. It contains a variety of active compounds, including eugenol, clove acids, and carvacrol, that give it potent antimicrobial properties.

**Procedure to use Material:** Insta Relief is a beneficial and potent natural remedy, which has been used to treat a variety of oral issues. It has antibacterial and antifungal properties that make it effective for dental pain, such as toothache. It may also help to reduce gum inflammation and promote healing. The active compounds in clove oil, eugenol and clove acids, are capable of killing a variety of oral bacteria. This makes it effective for treating oral infections, such as gum disease. It is usually applied directly to the gums for pain relief. The oil will quickly be absorbed by the gums and provide pain relief. You can also use clove oil as a mouthwash or toothache drops. It can also be used as a mouthwash to treat oral pain. Mix a few drops of clove oil with water and use it as you would a regular mouthwash. You can also use it as a toothache drop. Mix a few drops of it with a few drops of water and place the mixture in a clean dropper. Use the dropper to apply the clove oil to the gum line and teeth.

**Storage :** Keep Away from Direct Sunlight, Keep Away from Moisture, Always keep the Bottles and packets tightly closed once opened. Store at Temp. between 5°C to 30°C Shelf life not affected after opening if stored in a cool dry place, If not to be used for a longer duration, please keep them in refrigerator Allow material to reach room temperature prior to use

**Warning:** Avoid prolonged or repeated contact with skin (allergic contact dermatitis), and eyes. Avoid prolonged inhalation. Do not take internally. Eye contact: if the product comes in contact with eyes it may cause irritation. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention. Skin Contact: if comes in contact with Skin, it may cause skin sensitization, If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. It should not be used with patients who have a history of severe allergic reaction to any of the components

**Presentation:** 1 x 15 ml RefNo. MD02-01

**Shelf Life:** 3 years from date of manufacturing.

**Disposal Conditions:** To be disposed off based on the prevailing rules of the country in which it is used.

**NOTE: TO BE USED BY DENTAL PROFESSIONAL ONLY**

Manufactured By:  
**Golden Bird International**  
Plot No. A-4 Vishwakarma  
Industrial Estate, Village Panchli  
Khurd Meerut-250002 India, info@indiagbi.com  
www.dentakris.com, GOLDEN-Meeru-UP/M/MD/005962  
®-Registered Trademark Only

Symbol for "BATCH CODE" Symbol for "USE BY" Symbol for "NON STERILE" Symbol for "KEEP DRY" Symbol for "CAUTION"  
 SYMBOL FOR EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE Symbol for "KEEP AWAY FROM SUNLIGHT"  
 Symbol for "KEEP OUT OF THE REACH OF CHILDREN" Symbol for "MANUFACTURER NAME & ADDRESS" Symbol for "CONSULT INSTRUCTION FOR USE" Symbol for "CATALOGUE NUMBER" Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"  
 Symbol for "TEMPERATURE LIMITATION"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.